

Product introduction

Dental Photographic Mirror & Contrastors



ORTAM ORTHODONTICS

Photographic Mirror Features

- Durable thick material with rounded edges for patient comfort
- Thick material lets clinicians to expand cheeks by the rigid mirror while taking pictures.
- Not coated. One-piece solid Stainless Steel raw material polished by skilled technicians.
- Sterilizable by all methods.
- Allows taking good quality pictures with the conventional digital cameras.

Photographic Mirror Code #'s

Photographic Mirrors

Adult Occlusal/Child Occlusal

801-020



25 \$

Buccal Lateral

801-021



Adult Occlusal/Child Occlusal

801-022



Adult Buccal/Adult Occlusal

801-023



Photographic Mirror Code #'s


Lingual Upper/Lingual Lower

801-024



Child Buccal/Adult Buccal


801-025



25 \$

Child Occlusal/Child Buccal


801-026



Photographic Mirror Recommendations

- Before sterilization wipe and dry with soft cotton fabric to avoid spots
- Do not thumble mirrors with other instruments after use, this will cause scratches on the surface
- After long-term use, wiping dental stone-alcohol mixture via cotton piece can solve minor scratches and staining.
- Non-intended scratches and stains remaining on the mirrors are under guarantee of the manufacturer; they can be re-polished at 1/2 price for the lifetime. (Guarantee excludes freight charges and deep scratches)
- To avoid fogging while taking pictures, simply apply heat to back surface of the mirror. (Caution: Do not over heat to avoid cheek impingement.)

Intermaxillary Lip Bumper



Lateral Views taken By SS Mirror

Patient Treated By Spring Jumper Lateral Views Taken By SS Mirror

Pretreatment November 2007



Post Treatment May 2010

Conventional CL II Therapy Lateral Views taken By SS Mirror


Pre Treatment November 2007




Post Treatment November 2010

Dental Photo Contrastors

Occlusal Contrastors



25 \$



15 \$

Lip Depressor Contrastor Std.

25 \$

5545




7 \$



Plexiglass Standart
Lip depressor
45-55mm

Lip Depressor Contrastors

5550



25 \$

6055




Photo Retractors

Auto-clavable, up to 130°C



7 \$

G311-11

1 pc/pk

Lateral



12 \$

G311-12

2 pc/pk

Occlusion-Large



9 \$


G311-13

2 pc/pk

Occlusion-Small



ORTAM ORTHODONTICS



E-mail : catagun@yahoo.com

Dr. Çınar Atagün

TEL : +90 232 4222871

FAX : +90 232 4638195

Contact Information